

VisuMax Defining the pulse rate in refractive surgery



We make it visible.



// VisuMax MADE BY CARL ZEISS

Its remarkable features are its precision and innovative detail.

New trends in modern corneal surgery.

With the VisuMax[®], Carl Zeiss is significantly shaping the world of refractive surgery. This ground-breaking laser system employs high-performance femtosecond laser technology and is characterized by its outstanding cutting precision, unsurpassed speed and gentle treatment technique. The VisuMax is thus the ideal platform for therapeutic and refractive applications of cutting-edge corneal surgery, including Flap, Keratoplasty, Incision for ICR and ReLEx[®].

With ReLEx, VisuMax heralds a paradigm shift in refractive surgery: minimally invasive laser vision correction. It completes the unparalleled range of innovative surgical possibilities and creates the ideal preconditions for tapping into new patient groups.

VisuMax is the logical enhancement of the product range for refractive surgeons – and marks another step into the future of corneal surgery, the progress of which Carl Zeiss has been shaping for more than 20 years.

VisuMax applications

Precision in all facets

The VisuMax[®] is extraordinary in all major applications of state-of-the-art corneal surgery. The innovative femtosecond laser system brings together perfectly coordinated components to provide maximum cutting precision, efficiency, predictability and comfort.

ReLEx

With the minimally invasive refractive procedure, ReLEx[®] smile, VisuMax enables a new, flapless technique. In a single step, the femtosecond laser creates the refractive lenticule, and the access through which the lenticule is extracted. Without ablation and without a flap.

Flap

In Femto-LASIK and Laser Blended Vision for treating presbyopic patients, the VisuMax stands out as a flap cutter. It provides predictable flap thickness and adjustable geometries.

Keratoplasty

VisuMax offers a broad spectrum of corneal transplant procedures, including lamellar and penetrating keratoplasty. High-precision cutting quality and rapid incision progress enable the preparation of precision corneal grafts and ideal preparation of the recipient's cornea.

Incision for ICR

The VisuMax also impresses with its femtosecond laser technology for the implantation of intracorneal rings. It even permits inclined cutting geometries and ring tunnels smaller than 360° and offers a previously unachieved degree of flexibility when defining the tunnel parameters.



VisuMax strengths

Building blocks of state-of-the-art femtosecond technology



A contact glass as ingeniously designed as the cornea

The surface of the human cornea is curved. Thus, Carl Zeiss contact glasses are curved way, too. The three different sizes available (S, M, L) ensure an optimal fit to the anatomy of the eye. The cornea is not forced into a deformed planar, non-physiological shape – and artifacts are avoided in the cutting result, as is unnecessarily high IOP.



Maximum cutting precision

High-precision optics from Carl Zeiss provide an extremely focused laser beam. The result: minimum laser pulse energy at a high pulse frequency for hitherto unsurpassed incision control – at precisely the desired depth in the cornea, even with threedimensional, curved incisions.



Brilliant visual control

The integrated, high-quality ZEISS surgical microscope, including digital video camera for live recording of the surgical procedure ensures precise and complete control of each treatment step.



A smart unit

The ergonomically pivoting patient supporting system ensures maximum comfort. The patient's position is continuously monitored during treatment and the sturdy yet comfortable patient supporting system is automatically adjusted during surgery.

Achieving objectives with intuition

The VisuMax can be easily controlled via its touchscreen and intuitive software. An interactive wizard supports the surgeon during all steps.



Efficiency paying off

With a laser pulse frequency of 500 kHz, the VisuMax realizes short treatment times. This means more comfort for both physician and patient. In addition, the user benefits from a more efficient workflow and a higher throughput of satisfied patients.



Slit illumination for greater detail

As a universal workstation for corneal surgery, the VisuMax features integrated slit illumination for monitoring and control immediately after the respective treatment, without the patient having to switch places or change position.

ReLEx For a new flapless treatment

With ReLEx[®], the VisuMax[®] enables the fusion of cutting-edge femtosecond technology and precise lenticule extraction to provide minimally invasive vision correction. A refractive lenticule is created in the intact cornea and removed via a small incision. Without ablation. Without a flap. The treatment is thus a **flapless**, **all-femto** and **single-step** procedure.

Flapless

ReLEx is turning the world of refractive surgery on its head. In ReLEx smile, a small incision is sufficient to extract the lenticule. The minimally invasive incision also implies fewer transected nerves and significantly reduced incidence of dry-eye syndrome. The risk of infections, epithelial ingrowths or flap complications can be reduced. The smaller incisions enable better healing of the epithelium.

All-femto

With the VisuMax, ReLEx relies completely on femtosecond technology. The unique laser vision correction procedure creates the pre-calculated lenticule in the cornea with precision and predictability. Nomograms or fluence tests are not required, intraoperative ambient conditions or individual corneal characteristics have virtually no influence on the reproducibility of the lenticule. Physicians benefit from excellent predictability, particularly when correcting higher refraction values.

Single-step

With ReLEx, both the lenticule and the access incision are created in a single treatment step. Unlike Femto-LASIK, therefore, only one surgical procedure needs to be planned. There is no need for the patient to switch places. The result is more efficient workflows and shorter treatment times. For patients, the procedure is a much less stressful experience.



Outstanding results

- Small incision of less than 4 mm
- Side-cut length up to 80% shorter and
- Cap incision area up to 30 % smaller than for a Femto-LASIK flap
- Implies lower incidence of dry eye syndrome and less nerve transection due to smallest of incisions with no flap
- Less risk of infections, epithelial ingrowths and flap complications
- Reproducibility of the lenticule, irrespective of individual corneal characteristics and ambient conditions
- Excellent predictability, particularly when correcting higher refraction values
- Efficient treatment process without patient having to switch places

ReLEx smile

Small Incision Lenticule Extraction



In a single step, the VisuMax creates a refractive lenticule and a small incision of less than 4 mm in the otherwise intact cornea - almost irrespective of the ambient conditions or corneal structure.



The lenticule is removed through the small incision. The disruption to the biomechanics of the cornea is minimal. No flap is cut.



The removal of the lenticule alters the shape of the cornea, thus achieving the desired refractive change.

Flap For best possible conditions

For Femto-LASIK and the treatment with Laser Blended Vision, VisuMax[®] means above all one thing: high-precision flaps. Combined with the excimer laser MEL 80[™] and the treatment planning station CRS-Master[®], it offers an optimally coordinated system solution for refractive laser surgery – for a convenient workflow, efficient patient management and the best possible results.

The combination of precision and efficiency

- High-precision flaps due to highperformance femtosecond technology
- High reproducibility and consistency of flap thickness
- Easy repositioning of the flap
- Optimum workflow due to perfectly coordinated system components
- Smooth, finely structured surface of the stromal bed
- Prevents unnecessarily raised IOP and thus temporary losses of vision and trauma, due to the anatomically curved contact glass and the nonscleral suction.

Pivotable patient supporting system – Optimum workflows for greater comfort

The shared use of the pivotable patient supporting system by the VisuMax and the MEL 80 saves the patient from having to move from one treatment location to another. The patient experiences the surgery as an integrated process. Unnecessary waiting periods are minimized and treatment efficiency increased.



Femto-LASIK

Laser-in-situ-Keratomileusis



The VisuMax femtosecond laser creates the flap.



MEL 80 – All you need for optimum results

All the parameters of this high-precision excimer laser are geared to increasing efficiency, achieving optimum treatment results and rapid visual recovery. Key factors here are the extremely high ablation speed, customized treatment planning with the optional CRS-Master, the high-performance eyetracker system and eye registration.

CRS-Master – For a truly individual treatment

As a state-of-the-art treatment planning tool, the CRS-Master transfers data from the wavefront diagnosis and corneal topography into the MEL 80. This additional incorporation of patient-related data makes it possible to create a complete, individual eye profile for Femto-LASIK or Laser Blended Vision that is perfectly tailored to the patient.



The patient moves to the MEL 80 excimer laser.



The flap is manually opened and folded back to expose the deeper corneal layer (stroma) beneath.



The MEL 80 excimer laser ablates the pre-calculated corneal tissue point by point.



The upper corneal layer is then repositioned following the refractive correction.

Keratoplasty For high-precision tissue grafts

With the Keratoplasty option, the VisuMax[®] quickly converts into a state-of-the-art workstation for corneal grafts. The excellent cutting quality and laser control enable smooth lamellar and circular incision areas for high-precision results.

Broad spectrum for keratoplasty

The Keratoplasty option for the VisuMax comprises functions specially developed for corneal grafts and their optimization. High-precision and rapid incision progress, combined with a high level of reproducibility, mean that the VisuMax can be used for the three most important corneal graft procedures:

- Penetrating Keratoplasty (PKP)
- Deep Anterior Lamellar Keratoplasty (DALK)
- Descemet's Stripping Endothelial Keratoplasty (DSEK)

Perfect tissue grafts

The unique adapter attached to the headrest of the patient supporting system serves as an ideal work platform for preparing the corneal donor graft as well as the recipient's cornea.

Specially developed contact glass (type KP)

The curved contour of the contact glass prevents any unnecessary compression of the corneal tissue. It is also compatible with most artificial anterior chambers.



The practical adapter provides a robust and sterile work surface for preparing the corneal graft.



Contact glass (type KP) for preparing the donor cornea

The option "Endothelial Keratoplasty" is not intended for sale in the United States.



Precise results

- High-precision cutting quality in anterior lamellar and endothelial Keratoplasty
- Penetrating Keratoplasty with perfect fit for donor and recipient cornea
- Precisely predictable incision pattern for greater reliability when preparing thin grafts
- Small spot distance for excellent cutting quality and easy separation of the tissue
- Flexible adjustment of the cutting parameters on the VisuMax

Maximum efficiency for optimum workflows

- 500 kHz laser pulse frequency for faster, more precise treatments and shorter cutting times (typically less than 60 seconds), even for very deep cuts
- High-quality surgical microscope for all treatment phases
- Practical adapter for preparing the donor cornea
- Special contact glass (type KP) compatible with most artificial anterior chambers



View through the surgical microscope of the VisuMax: Separation of the cut lamellae from the recipient cornea

Incision for ICR

For flexible access

When implanting intracorneal ring (ICR) segments, surgeons benefit from the unique advantages of the VisuMax[®] femtosecond laser. The Incision for ICR option offers the possibility, for the first time, to create even inclined cutting geometries and tunnel segments between 90° and 270°. The repeatedly proven femtosecond laser technology ensures not only high-precision cutting quality but also previously unattainable degrees of freedom when defining the tunnel parameters.



Freely variable cutting parameters: for the first time, an incision can even be created parallel to the posterior surface of the cornea.

Degrees of freedom redefined

The VisuMax allows the corneal tunnels for the implantation of intracorneal ring (ICR) segments to be prepared quickly and easily, in the precise position and depth required. The wide range of adjustable parameters and the various combinations thereof give surgeons a unique degree of flexibility.

Tailor-made segments

For the first time, it is also possible to create tunnels with an arc angle of less than 360° using a femtosecond laser: With segments of between 90° and 270° partial tunnels can be designed individually, with precision and a high degree of flexibility. This means that individual intracorneal ring segments can also be implanted, without creating unnecessary incisions in the corneal tissue. Thus, tunnels can also be created for mixed ring segments with different geometries.



New treatment alternatives offer individual benefits

- The option to select between 0, 1 or 2 trapezoidal access incisions facilitates flexible tunnel access
- Width and inclination of the tunnel can be freely defined and precisely adjusted to the individual corneal shape and the applied ring geometry
- Seamlessly integrated into the user interface of the VisuMax, the ICR option provides maximum ease of use
- Rapid and intuitive entry of the necessary parameters

- The possibility to save user-defined cutting geometries as reusable templates increases workflow efficiency
- Reliability enhancement by graphic visualization of the parameter selection and automatic consistency check of the input parameters
- The unique contact glass is modelled on the natural shape of the corneal surface and thus contributes to patientfriendly treatment methods.
- Excellent control of tunnel preparation and ICR implant insertion as well as complete video documentation using the high-quality ZEISS surgical microscope

All common ICR products are supported.

The option Incision for ICR is not intended for sale in the United States.



Technical data

Installation and operating instructions

VisuMax femtosecond laser system

visuwax remitosecona laser system			
System components	Patient supporting s	system, including platform	
	Integrated uninterru	iptible power supply (UPS)	
	Surgical microscope	with additional slit illumination	
	Integrated with digital recording video camera		
Laser parameters	Wavelength	1043 nm	
	Pulse duration	220–580 fs	
	Laser pulse rate	500 kHz	

Installation and set-up conditions

Weight	870 kg (including patient supporting system, platform, UPS)	
Recommended space requirement	4.40 m x 3.80 m (standalone)	
Electrical connection	100–240 V, 50/60 Hz, max. 16 A	
	Separately fused circuit	

Operating conditions	
Room temperature	18 to 25 °C
Atmospheric humidity	30 to 70 %
Accessories	Single-use contact glasses Treatment Pack (sizes S / M / L and type KP) Keratoplasty adapter for patient supporting system



Your local contact:

Argentina Carl Zeiss Argentina S.A. Calle Nahuel Huapi 4015 / 25 C1430 BCO Buenos Aires Argentina Phone: +54 11 45 45 66 61 bruzzi@zeiss.com.ar

Australia Carl Zeiss Ptv Ltd

Tenancy Office 4, Level 1 40-52 Talavera Road North Ryde NSW 2113 Australia Phone: +61 2 9020 1333 med@zeiss.com

Austria Carl Zeiss GmbH Laxenburger Str. 2 1100 Vienna Austria Phone: +43 1 79 51 80 austria@zeiss.org

Belgium

Carl Zeiss NV-SA Ikaroslaan 49 1930 Zaventem Belgium Phone: +32 2 719 39 11 info@zeiss.be

Brazil

Carl Zeiss do Brasil Ltda. Av. Naçoes Unidas, 21711 CEP04795-100 São Paulo Brazil Phone: +55 11 5693 5521 medbrasil@zeiss.org

Canada

Carl Zeiss Canada Ltd. 45 Valleybrook Drive Toronto, ON M3B 2S6 Canada Phone: +1 800 387 8037 micro@zeiss.com

China

Carl Zeiss Shanghai Co. Ltd. 1/f., Ke Yuan Building 11 Ri Yin Nan Road Waigaoqiao Free Trade Zone 2005 Yang Gao Bei Road Shanghai 200131 China Phone: +86 21 5048 17 17 sro@zeiss.com.cn

CE₀₂₉₇

Manufacturer:

Carl Zeiss Meditec AG Goeschwitzer Strasse 51–52 07745 Jena Germany www.meditec.zeiss.com/VisuMax Czech Republic Carl Zeiss spol. s.r.o. Radlická 14/3201 150 00 Prague 5 Czech Republic Phone: +420 233 101 221 zeiss@zeiss.cz

France Carl Zeiss Meditec France SAS

60, route de Sartrouville 78230 Le Pecq France Phone: +33 1 34 80 21 00 med@zeiss.fr

Germany

Carl Zeiss Meditec VG mbH Carl-Zeiss-Strasse 22 73446 Oberkochen Germany Phone: +49 7364 20 6000 vertrieb@meditec.zeiss.com Surgical Ophthalmology: Phone: +49 800 470 50 30 iol.order@meditec.zeiss.com

Hong Kong

Carl Zeiss Far East Co. Ltd. Units 11-12. 25/F Tower 2, Ever Gain Plaza No. 88 Container Port Road Kwai Chung Hong Kong Phone: +852 2332 0402 czfe@zeiss.com.hk

India

Carl Zeiss India Pvt. Ltd. 22. Kensington Road Ulsoor Bangalore 560 008 India Phone: +91 80 2557 88 88 info@zeiss.co.in

Italy Carl Zeiss S.p.A. Viale delle Industrie 20

20020 Arese (Milan) Italy Phone: +39 02 93773 1 post@zeiss.it

Japan

Carl Zeiss Meditec Japan Co. Ltd. Shinjuku Ku Tokyo 160-0003 22 Honchio-Cho Japan Ophthalmic instruments: Phone: +81 3 33 55 0331 medsales@zeiss.co.jp Surgical instruments: Phone: +81 3 33 55 0341 cmskoho@zeiss.co.jp

Malaysia

Carl Żeiss Sdn Bhd. Lot2, Jalan 243/51 A 46100 Petaling Jaya Selangor Darul Ehsan Malaysia Phone: +60 3 7877 50 58 malaysia@zeiss.com.sg

Mexico

Carl Zeiss de México S.A. de C.V. Avenida Miguel Angel de Quevedo 496 04010 Mexico City Mexico Phone: +52 55 59 99 0200 cz-mexico@zeiss.org

Netherlands

Carl Zeiss B.V. Trapezium 300 Postbus 310 3364 DL Sliedrecht Netherlands Phone: +31 184 43 34 00 info@zeiss.nl

New Zealand

Carl Zeiss NZ Ltd 15B Paramount Drive P.O. Box 121 - 1001 Henderson, Auckland 0650 Ph: +64 9 838 5626 med@zeiss.com

Poland Carl Zeiss sr

Carl Zeiss sp. Z o.o. ul. Lopuszanska 32 02-220 Warsaw Poland Phone: +48 22 858 2343 medycyna@zeiss.pl

Singapore Carl Zeiss Ptd. Ltd. 50 Kaki Bukit Place Singapore 415926 Singapore Phone: +65 6741 9600 info@zeiss.com.sq

South Africa Carl Zeiss (Pty.) Ltd. 363 Oak Avenue Ferndale Randburg 2194 South Africa Phone: +27 11 886 9510 info@zeiss.co.za South Korea Carl Zeiss Co. Ltd. Seoul 121-828 Mapo-gu 141-1, Sangsu-dong 2F, BR Elitel Bldg. South Korea Phone: +82 2 3140 2600 korea@zeiss.co.kr

Spain

Carl Zeiss Meditec Iberia S.A. Ronda de Poniente, 15 Tres Cantos 28760 Madrid Spain Phone: +34 91 203 37 00 info@zeiss.es

Sweden

Carl Zeiss AB Tegeluddsvaegen 76 10254 Stockholm Sweden Phone: +46 84 59 25 00 info@zeiss.se

Switzerland Carl Zeiss AG Feldbachstrasse 81 8714 Feldbach Switzerland Phone: +41 55 254 7534 med@zeiss.ch

Thailand Carl Zeiss Thailand Floor 8, Thosapol Land Building 2 230 Ratchadapisek Road Huaykwang, Bangkok 10310 Thailand Phone: +66 2 2 74 06 43 thailand@zeiss.com.sg

United Kingdom Carl Zeiss Ltd. 15-20 Woodfield Road Welwyn Garden City Hertfordshire, AL7 1JQ United Kingdom Phone: +44 1707 871200 info@zeiss.co.uk

United States of America Carl Zeiss Meditec, Inc. 5160 Hacienda Drive Dublin, CA 94568 USA Phone: +1 925 557 4100 info@meditec.zeiss.com

from the current status of approval of the product in your country. Please contact our regional representative for more information. The contents of the brochure may differ from the current status of approval of the product in your country. Please contact our regional representative for Subject to change in design and scope of delivery and as a result of ongoing technical development. Printed on elemental chlorine-free bleached paper. © 2012 by Carl Zeiss Meditec AG. All copyrights reserved. CZ-VI/2012 Printed in Germany Publication No: 000000-1980-341

г

L

٦